

ASCENDING AORTIC PROSTHESES

CARBOMEDICS CARBO-SEALTM & CARBO-SEAL VALSALVATM







CARBOMEDICS CARBO-SEAL[™] & CARBO-SEAL VALSALVA[™] Over 30 years of experience in mechanical heart valves

Carbomedics Carbo-Seal AAP includes the Carbomedics Standard Aortic Valve, a fully rotatable valve with a very low rate of thromboembolic events, favorable hemodynamics, and negligible postoperative structural failures.¹⁻⁶



FLEXIBILITY, HANDLING AND NATURAL SHAPE TO COEXIST IN A BETTER WAY WITH NATURE

AAP: Aortic Ascending Prostheses



Flexibility and Handling

- Thin, pliable, woven polyester from Vascutek[®] requires no pre-clotting.
- Easier handling and suturing in comparison to bulkier velour materials.¹²
- Gelatin coated and crimped structure for optimal flexibility.
- Ultra-low porosity fabric and gelatin sealing result in less leakage, weeping and blushing.
- Resists fraying and minimizes suture hole bleeding.¹⁰
- Orientation markers provide easy visual suture positioning.





ATTIONS



CARBOMEDICS CARBO-SEAL VALSALVA[™] A better way to coexist with nature

SINUS OF VALSALVA

- Vertical orientation of pleats facilitates coronary anastomosis.7
- Reproduces the native sinus, reducing stress on the coronary anastomosis.⁷
- Sinus design encourages natural formation of systolic vortex.89
- Low profile taper to reduce stress on the coronary ostia buttons.⁷
- Sealed graft with no pre-clotting requirements.
- Full-sized standard aortic valve provides favorable hemodynamics.^{1,2}



Natural blood flow



Gelatin promotes natural healing¹⁰

Carbomedics Carbo-Seal AAP's Gelweave[™] graft is infused with minimally crosslinked gelatin that does not alter the healing process, encouraging a secure neo-intimal attachment with reduced inflammatory response.¹⁰

100 80

60

40

20

0



The histology at three months shows a typical, well-incorporated Dacron graft. The pseudointima is thin and is not proliferating.



Unlike collagen-coated grafts which are enzymatically metabolized, the Carbomedics Carbo-Seal's gelatin sealant is biodegraded by hydrolysis and rapidly absorbed within 14 days.*

* CORCYM data on file

User friendly instrumentation

Carbomedics Carbo-Seal and Carbomedics Carbo-Seal Valsalva AAP come on a holder assembly to facilitate initial proximal placement.



Disposable cautery included inside every implant packaging.



Dual-ended.

	100	-2-	Convers 1	
	1000	-8	Calcula I	 and a
	-	-8	Concerned a loss	
	100	*	(before)	 R
	(1.5 million	
tester.	-	<u> </u>		

CARBOMEDICS **CARBO-SEAL**

ASCENDING AORTIC PROSTHESIS (AAP) Sizes 21-33 mm

PRODUCT ORDERING INFORMATION

Valve Size	Valve Orifice Area (cm²)	Graft Inner Diameter (mm)	Minimum Graft Lenght (cm)	Catalog Number
21	2.07	24	10	AP-021
23	2.56	26	10	AP-023
25	3.16	28	10	AP-025
27	3.84	30	10	AP-027
29	4.44	32	10	AP-029
31	4.44	34	10	AP-031
33	4.44	34	10	AP-033



CARBOMEDICS **CARBO-SEAL VALSALVA**

ASCENDING AORTIC PROSTHESIS (AAP) Sizes 21-29 mm

PRODUCT ORDERING INFORMATION

Valve Size	Valve Orifice Area (cm²)	Graft Inner Diameter (mm)	Maximum Sinus Inner Diameter (mm)	Sinus Region Lenght* (mm)	Minimum Graft Lenght (cm)	Catalog Number
21	2.07	24	32	24	10	CP-021
23	2.56	26	34	26	10	CP-023
25	3.16	28	36	28	10	CP-025
27	3.84	30	38	30	10	CP-027
29	4.44	32	40	32	10	CP-029

* Based on cadaveric anatomic studies, it was determined that the length of the sinus region should equal the bore diameter of the straight portion of the graft.

ACCESSORIES ORDERING INFORMATION

CODE	ARTICLE	DESCRIPTION	
TR-101	Empty tray	1 empty tray	100 A
VS-200	Sizer set	4 sizers 19-21mm, 23-25mm 27-29mm, 31-33mm	
AR-150	Rotators set	Rotators set	
VT-100	Occluder tester	10 disposable occluder tester (provided sterile)	

REFERENCES

- 1. Bouchard et al., Twenty-Year Experience With the CarboMedics Mechanical 5. Ryuji Tominaga et al., A 10-year experience with the Carbomedics cardiac 10. Drury et al., Experimental and Clinical Experience with a Gelatin Valve Prosthesis. Ann Thorac Surg 2014;97:816-23
- 2. Aagard et al., Midterm Evaluation of Hemodynamics of the Top Hat 6. Carrier, Fifteen Years of Clinical and Echocardiographic Follow Up with the 11. Cavallaro A. et al., 1986. Pretreatment with gelatin of Dacron Grafts. Supraannular Aortic Valve. Asian Cardiovasc Thorac Ann 2010;18:1–5
- 3. Langley et al., Replacement of the proximal aorta and aortic valve using a 7. composite bileaflet prosthesis and gelatin-impregnated polyester graft (Carbo-seal): early results in 143 patients, SM. J Thorac Cardiovasc Surg 1999:118:1014-20
- 4. De Paulis R. et al., Opening and closing characteristics of the aortic valve after valve-sparing procedure usign a new aortic root conduit, Ann Thorac 9. Surg 2001;72:487-494
- prosthesis. Ann Thorac Surg 2005;79:784-9
- CarboMedics Heart Valve. J Heart V Dis 2006 Jan; 15(1): 67-72
- De Paulis et al., One-year appraisal of a new aortic root conduit with sinuses of Valsalva, J Thorac Cardiovasc Surg 2002 Jan;123(1):33-9
- 8 De Paulis et al., Analysis of Valve Motion After the Reimplantation Type of Valve-Sparing Procedure (David I) With a New Aortic Root Conduit, Ann Thorac Surg 2002;74:53-7
 - De Paulis et al., A New Aortic Dacron Conduit for Surgical Treatment of Aortic Root Pathology, Ital Heart J 2000, 1 (7),457-463
- Impregnated Dacron Prosthesis, Ann Vasc Surg, 1987, 7,542-547
- Experimental Research. Policlinico Sez Chir. Vol. 93 pp. 275-283 [article in italian]
- 12. Kadoba et al., Experimental comparison of albumin-sealed and gelatin-sealed knitted Dacron conduits, J Thorac and Cardiovasc Surg 1992;103:1059-67

INTENDED USE/INDICATIONS

The Carbomedics Ascending Aortic Prostheses are intended for use in open heart surgery for simultaneous replacement of the ascending aorta and the aortic valve in cases of aneurysm, dissection, or other disease conditions of the aorta combined with disease or degeneration of the aortic valve.

KEY CONTRAINDICATIONS

The Carbomedics Ascending Aortic Prostheses are contraindicated in patients unable to tolerate long term anticoagulation therapy or for whom this type of therapy is difficult to maintain.

KEY WARNINGS

For single use only. Do not attempt to clean, resterilize, or reuse any prosthesis that has been in contact with organic blood or tissue. This product should not be implanted in patients who exhibit sensitivity to polyester or material of bovine origin. Do not preclot.

TOP POTENTIAL SIDE EFFECTS

The risks or potential adverse events associated with the use of prosthetic aortic heart valves include, but Adverse events potentially associated with the use of prosthetic aortic heart valves include: angina; cardiac

arrhythmia; endocarditis; heart failure; hemolysis; hemolytic anemia; hemorrhage; myocardial infarction; prosthesis nonstructural dysfunction; prosthesis perivalvular leak; prosthesis regurgitation; prosthesis structural dysfunction; prosthesis thrombus; stroke; thromboembolism. Adverse events potentially associated with the use of vascular grafts include: aneurysm; clinical reactions to collagen implantation have been described as infrequent, mild, localized, and self limiting; embolism; hemorrhage; infection; occlusion; pseudoaneurysm; seroma. It is possible that these complications could lead to: reoperation, explant, permanent disability, death.

MRI conditional

For professional use. Instructions for Use are available upon request through the manufacturer's website. Not approved in all geographies. Consult your labeling



Manufactured by:

Corcym S.r.l. Via per Crescentino sn 13040 Saluggia (VC) Italy Tel: +39 0161 487800

Corcym S.r.l. previously Sorin Group Italia S.r.l.



corcym .com

© 2021 all rights reserved

CC-MK-00004 B