

CARDIOXYL® CARDIOXYL


Tresse de polyester enduite de silicone

Chirurgie cardiopédiatrique

Pr Kreitmann

Hôpital de la Timone, Marseille, France

A choisi Cardioxyl® pour la fermeture
des CIV* à points séparés



- Pledgets et aiguilles adaptés à la taille des patients
- Confort d'utilisation grâce à une souplesse idéale et une absence de mémoire de forme
- Fil tressé assurant une solidité remarquable et une excellente tenue au nœud
- En 2 couleurs (vert et blanc) pour une organisation facilitée des sutures

Référence	Couleur	USP	Déc	Nombre d'aiguille	Courbure aiguille	Type aiguille	Taille aiguille	Fil	Pledget	Conditionnement
73P10B	3 vert +	5/0	1	Double aiguillée	1/2	Ronde	10 mm	60 cm	Weavenit 2 x 3 mm	Boîte de 12 unités
73P10A						KL	13 mm			
73P07A	3 blanc	6/0	0.7			Ronde	10 mm			

*CIV : Communication Interventriculaire



CARDIOXYL® CARDIOXYL


Silicone coated polyester braid

Paediatric Cardiac surgery

Pr Kreitmann

Timone Children's Hospital, Marseille, France

Has chosen Cardioxyl® for VSD*
closure with interrupted sutures



- Pledgets and needles adapted to the size of the patients
- Comfortable use thanks to increased flexibility and a lack of thread memory
- Braid allowing exceptional strength and knot holding
- Available in 2 colors (green and white) to facilitate sutures organisation

Reference	Color	USP	Metric	Needle number	Needle curvature	Needle point	Needle size	Thread	Pledget	Packaging
73P10B	3 green + 3 white	5/0	1	Double armed	1/2	Taper	10 mm	60 cm	Weavenit 2 x 3 mm	Box of 12 units
73P10A						Tapercutting	13 mm			
73P07A		6/0	0.7			Taper	10 mm			

*VSD: Ventricular Septal Defect

