

RECTAL CANNULA



VECTEC®
SINGLE USE
LAPAROSCOPIC
DEVICES



Single use cannula

- Allows excellent recto sigmoid exposure.
- Facilitates double safety testing.
 - 1- Air test
 - 2- Liquid dye test



Reference	Designation	Shaft diameter	Tip diameter	Tip length	No. unit per box
RCC 320	Canule rectale	4 mm	32 mm	50 mm	10

Product description :

- A distal cannula with:
 - 1- an ergonomic and atraumatic tip for safe progression into the rectal and sigmoid cavities.
 - 2- a flat rear part to the tip good for safety tests (air/liquid).
- A small-diameter stainless steel shaft which is flexible to adapt to the patient's anatomy.
- An ergonomic handpiece.
- T tubing:
 - 1- a length of small diameter silicone tubing with a stop-cock valve and luerlock for air/liquid injection, to be located within the sterile zone.
 - 2- a length of large diameter tubing with a red clamp for waste.



MD Class IIa

See instructions for use

CE
0459

G
MED NF EN ISO 13485



Péters
SURGICAL

Manufactured by:

Péters Surgical
42 rue Benoit Frachon
93013 Bobigny Cedex - FRANCE
Tel: +33 (0)1 48 10 62 62
Fax: +33 (0)1 48 91 22 99
commercial@peters-surgical.com
www.peters-surgical.com